

```

1 'Roni Hogri 4.7.22
2 'Script to automatically analyze calcium signals
3 '
4
5
6 'Imported data should be non-normalized (raw data from imageJ as txt file) - maximum of
  15 channels (cells) per file!
7 'f0 is defined locally from the last N (normally 30) subthreshold bins (1 s each).
8 'Filtered channels (DC and/or IIR) are marked by "*".
9 'If DCTC (DC time constant) or IIR frequency and/or are set to zero, DC remove and/or IIR
  filtration are skipped, respectively.
10 'After creating the smr file with the dF/F values, the process of looking for Ca2+ events
  begins:
11 'A threshold is calculated for each channel, and the times of suprathreshold events are
  marked in a level channel (1 per waveform channel).
12 'Threshold = baseline mean + multiplier x SD of baseline. Search for suprathreshold
  events starts with treatment onset (end of baseline)
13
14
15
16 'Outputs:
17 '1. Text file with normalized calcium signal intensities (dF/F values), after filtration
  (if set)
18 '2. Text file with values of all parameters measured + analysis info
19
20
21 '*****
22
23
24 main();
25
26 proc main();
27
28
29 '*****'User-defined variables:
30
31 var BaseStart := 0; 'when should the baseline period start (s)?
32 var BaseEnd% := 30; 'when should the baseline period end (s)?
33 var NumFBins% := 30; 'number of bins used for calculation of df/f
34 var HPFreq := 0.001; 'corner frequency to high pass with the IIR filter - higher value is
  stronger filter
35 var FilType% := 1; 'type of filter: 0 = butterworth, 1 = bessell
36 var DCTC := 0; 'time constant for DC remove process (s) - lower value is stronger filter.
  Set to 0 if you don't want to use it
37 var Mult := 2; 'multiplier of SD for threshold calculation
38 var MinDur% := 5; 'minimal duration of event (in whole s), i.e. all points during this
  period must be suprathreshold for an event to be detected
39 '(does NOT affect dF/F calculation!)
40
41 'suffix of file names for the report files created by program:
42 var FileEnd1$ := "drift_corrected.txt"; 'for file with all data points for all cells
  after drift correction (dF/F and filtration)
43 var FileEnd2$ := "SUMMARY.csv"; 'for file holding extracted results (peaks, areas, etc...)
44
45 'what to do when script finishes:
46 var ShowLast% := 5; 'which file to show at the end: 0 = smr, 1 = log, 2 = script, 3 = txt
  file; window must not be closed by Closewins%!
47 var CloseWins% := 1; 'set to 1 if you want to close all files except the script once
  program finishes
48
49
50
51 '*****
52 'END OF USER-DEFINED VARIABLES
53 '*****
54
55
56
57 '*****' General variables and processes
58
59 var vs% := App(3); 'view handle for the script
60 var vh%; 'view handle for smr file
61 var vl% := LogHandle(); 'handle for log file (for printing output)
62 var tl%; 'view handle for first txt file (value per s after processing)
63 var OGFNE$, OGFN$; 'full name of original txt file (including path), with and without
  file extension

```

```

64 var NSMRN$; 'fullname of new smr file with running F values
65 var NumWav%; 'number of valid waveform chans
66 var mt%; 'number of seconds (lines) in data file (equivalent to maxtime)
67 var FT$; 'for filter types
68 docase
69 case FilType% = 0 then FT$ := "Butterworth";
70 case FilType% = 1 then FT$ := "Bessel";
71 endcase
72
73
74 FrontView(vl%); EditSelectAll(); EditClear(); 'clear log
75
76 DefTxt(OGFN$, OGFNE$, NSMRN$, mt%, NumWav%, vs%); 'read txt file with values from imageJ,
count rows and columns
77
78 var OGVals[mt%][NumWav%]; 'for storing the raw values per s from txt file created by imageJ
79 var BaseMean[NumWav%]; 'mean value of baseline period per chan
80 var BaseSD[NumWav%]; 'SD value of baseline period per chan
81 var Thresh[NumWav%]; 'threshold value per chan
82 var EventChan[NumWav%]; 'identifier for marker channels (threshold crossing)
83 var EventNum[NumWav%]; 'number of suprathreshold events per chan
84 var MaxBins% := (mt% - BaseEnd%) / MinDur%; 'maximal possible number of suprathreshold
events per channel
85 var RiseTime[NumWav%][MaxBins%]; 'times of suprathreshold events - start
86 var FallTime[NumWav%][MaxBins%]; 'times of suprathreshold events - end
87 var EventDur[NumWav%][mt% - BaseEnd%]; 'event durations
88 var Peak[NumWav%][MaxBins%]; 'peak amplitudes of suprathreshold events
89 var AOT[NumWav%][MaxBins%]; 'area over threshold per event
90 var dFF0[mt%][NumWav%]; 'for storing the calculated dF/F values throughout the session
91
92 ImpText(OGFNE$, OGVals[[]], mt%); 'import raw values from txt file
93
94 CalF(NumFBins%, mt%, NumWav%, Mult, OGVals[[]], dFF0[[]]); 'calculate the df/f with a
running window
95
96 WriteNewFs(mt%, NumWav%, NSMRN$, dFF0[[]], vh%); 'write df/f values onto new smr file
97
98 if DCTC > 0 or HPFreq > 0 then 'if filter(s) set by user
99     ApplyFilters(NumWav%, HPFreq, FilType%, DCTC); 'filter out slow drifts in all wave
chans
100 endif
101
102 ThreshCalc(BaseMean[], BaseSD[], Thresh[], NumWav%, Mult, BaseStart, BaseEnd%);
'Calculate baseline mean, baseline SD, and threshold per chan
103
104 ThreshMark(mt%, Thresh[], NumWav%, MinDur%, EventChan[], BaseEnd%); 'Create event chan,
mark suprathreshold events in each equivalent wavchan
105
106 EventParams(RiseTime[[]], FallTime[[]], EventDur[[]], Peak[[]], AOT[[]], NumWav%,
EventChan[], mt%, MinDur%, Thresh[], BaseEnd%);
107 'For each event, find peak amplitude and area over threshold
108
109 'print results to 2 separate files:
110 PrintResults_dFF(mt%, vh%, NumWav%, FileEnd1$, t1%, CloseWins%, OGFN$); '1. txt file with
the corrected values per chan per s, no titles
111 PrintResults_Summary(vh%, HPFreq, Mult, NumWav%, MaxBins%, EventDur[[]], Peak[[]],
AOT[[]], Thresh[], MinDur%, FileEnd2$, DCTC, CloseWins%,
112 OGFN$, FT$); '2. csv file with peak, area, threshold, info
113
114
115 if CloseWins% = 1 then
116     FileClose(-1, 0);
117 endif
118
119 docase 'which window to show when program finishes running
120 case ShowLast% = 0 then FrontView(vh%); 'show smr data file
121 case ShowLast% = 1 then FrontView(vl%); 'show log
122 case ShowLast% = 2 then FrontView(vs%); 'show script
123 case ShowLast% = 3 then FrontView(t1%); 'show text file (all values after filtration)
124 endcase
125
126
127 end 'end of main procedure
128
129
130

```

```

131 'Start of sub-procedures
132 '-----
133
134
135 '-----'read txt file with values from imageJ, count rows and columns
136 proc DefTxt(&OGFN$, &OGFNE$, &NSMRN$, &mt%, &NumWav%, vs%);
137
138
139 var vot%; 'view of txt file with original pixel values
140 vot% := FileOpen("",8,8,"Please select a .txt file to import raw data from");
141
142 if vot% < 0 then 'if no file selected
143     FrontView(vs%); 'return to script and abort
144     halt
145 endif
146
147 OGFN$ := FileName$(1) + FileName$(2) + FileName$(3) + FileName$(4); 'full path of
148 original file, without file extention
149
149 NSMRN$ := OGFN$ + ".smr"; 'name of new smr data file with dF/F values
150 OGFNE$ := OGFN$ + ".txt"; 'original file, with file extention
151
152 var col[100]; 'original array for holding column (channel) values from first row (time bin)
153
154 var c%;
155 ReadSetup("", "");
156 c% := Read(col[]); 'fill array with values of first row, return number of
157 cells/columns/channels
158
158 if c% > 15 then 'Avoid exceeding Spike2 channel number limitation
159     Message("Too many channels (columns) in txt file!! Maximum of 15 channels allowed!");
160     halt
161 endif
162
163 resize col[c%]; 'resize the columns array according to the actual number of channels
164
165 var i%; 'loop index
166
167 for i%:=0 to c%-1 do
168     if col[i%] > 0 then
169         NumWav% += 1; 'Spike2 channel number is column number + 1
170     endif
171 next
172
173 var line$; 'for reading lines from txt file
174 ReadSetup("", "");
175 while Read(line$) >= 0 do 'keep going until line is empty
176     mt% += 1; 'count number of lines (time bins) in text file
177 wend
178 mt% += 1; 'add 1 to line count since the reading of the first line for column counting
179 skipped a line
180
180 FileClose(); 'close txt file so that next read starts from top
181
182
183 end
184
185
186
187
188 '-----'import raw values from txt file
189 proc ImpText(OGFNE$, &OGVals[][], mt%);
190
191
192 FileOpen(OGFNE$, 8); 'open original text file and read it
193 ReadSetup("\n", "\t");
194
195 var ii%; 'loop index
196
197 for ii%:=0 to mt%-1 do
198     Read(OGVals[ii%][]); 'in each line (time bin), store column (channel) values in
199 subarray
200 next
201
202 end

```

```

203
204
205
206 '::::::::::::::::::::'calculate the df/f with a running window
207 proc CalF(NumFBins%, mt%, NumWav%, Mult, OGVals[][], &dFF0[][]);
208
209 var FsArray[NumFBins%][NumWav%]; 'holds current F values in running window
210 var locmean, locSD, locthresh; 'local values that change as the window advances
211
212 var i%, ii%, jj%; 'loop indices
213
214 for i%:=0 to NumWav%-1 do 'for each cell
215     for ii%:=0 to NumFBins%-1 do 'for each time bin
216         FsArray[ii%][i%] := OGVals[ii%][i%]; 'the first window contains all the F values
of the first N bins
217     next
218 next
219
220
221 for i%:=0 to NumWav%-1 do 'for each cell
222
223     ArrSum(FsArray[][i%], locmean, locSD); 'calculate the "local" mean and SD of F values
in the current window
224
225     for ii%:=0 to NumFBins%-1 do 'calculate the dF/F of the first window
226         dFF0[ii%][i%] := (OGVals[ii%][i%] - locmean) / locmean;
227     next
228     locthresh := locmean + locSD * Mult; 'original threshold based on baseline period
229
230     for ii%:=NumFBins% to mt%-1 do 'for all remaining time bins
231         if OGVals[ii%][i%] < locthresh then
232             'if the current F value is lower than the threshold, add the new F value to the
window and move the window by 1 bin
233             for jj%:=0 to NumFBins%-2 do
234                 FsArray[jj%][i%] := FsArray[jj%+1][i%];
235             next
236             FsArray[NumFBins%-1][i%] := OGVals[ii%][i%];
237         endif
238         ArrSum(FsArray[][i%], locmean, locSD); 'calculate the mean and SD of F values in
the current window
239         dFF0[ii%][i%] := (OGVals[ii%][i%] - locmean) / locmean; 'calculate dF/F value for
this time point in this cell
240         locthresh := locmean + locSD * Mult; 'updated threshold for inclusion of next
point in the running window
241     next
242
243 next
244
245
246 end
247
248
249
250 '::::::::::::::::::::'import data from txt file into smr file
251 proc WriteNewFs(mt%, NumWav%, NSMRN$, dFF0[][], &vh%);
252
253 vh% := FileNew(7, 0, 1000, 1, mt%); 'create and name the new smr file
254 FileSaveAs(NSMRN$, -1, 0);
255
256 var ChanName$; 'for naming the created channels
257 var mch%; 'memory channel
258
259 var i%, ii%; 'loop indices
260
261 for i%:=0 to NumWav%-1 do 'for each imported channel
262     mch% := MemChan(1, 0, 1, 0); 'create temporary memory channel (waveform), with offset
of 0 V
263     ChanOffset(mch%, 0);
264     for ii%:=0 to mt%-1 do 'for each time bin in this channel
265         MemSetItem(mch%, 0, ii%, dFF0[ii%][i%]); 'Add waveform point to mem channel at
time ii%
266     next
267
268     ChanSave(mch%, i%+1); ChanShow(i%+1); 'save memory channel; show and name permanent
channel
269     ChanName$ := "Chan " + Str$(i%+1); ChanTitle$(i%+1, ChanName$);

```

```

270 ChanDelete(mch%); 'delete memory channel
271 next
272
273
274 Draw(0, MaxTime()); 'optimise axes
275 Optimise();
276
277
278 end
279
280
281
282 ''''''''''''''filter out slow drifts in all wave chans
283 proc ApplyFilters(NumWav%, HPFreq, FilType%, DCTC);
284
285 var ChanName$; 'filtration is indicated in channel name
286
287 var chl%; 'memory buffer identifier
288 var i%; 'loop index
289
290 for i%:=1 to NumWav% do 'for all waveform channels
291     ChanName$ := ChanTitle$(i%); 'name of original channel
292
293     if HPFreq > 0 then 'apply IIR filter?
294         chl% := MemChan(0, i%); 'new memory channel with same properties as original chan
295         IIRCreate(-1, 1, FilType%, 2, HPFreq); 'IIR High pass 0.001 Hz, second order
296         IIRApply(-1, chl%, i%, 0, MaxTime(), HPFreq); 'Apply created filter - IIR
butterworth order of 2, time constant 0.001 s
297         ChanDelete(i%); 'delete original channel
298         ChanSave(chl%, i%); 'save memory channel instead of original channel
299         ChanShow(i%); 'show new channel
300         ChanDelete(chl%); 'delete memory channel
301     endif
302
303
304     if DCTC > 0 then
305         ChanProcessAdd(i%, 2, DCTC); 'Use the 'DC remove' process (additional filtration)
306     endif
307
308     ChanName$ += "*"; ChanTitle$(i%, ChanName$); 'add * to channel name to indicate
filtration
309     Optimise(); 'optimize y axis to filtered signal
310
311 next
312
313
314 end
315
316
317
318 ''''''''''''''Calculate baseline mean, baseline SD, and threshold per channel
319 proc ThreshCalc(&BaseMean[], &BaseSD[], &Thresh[], NumWav%, Mult, BaseStart, BaseEnd%);
320
321 var i%; 'loop index
322
323 for i%:=0 to NumWav%-1 do 'for each cell
324     BaseMean[i%] := ChanMeasure(i%+1, 2, BaseStart, BaseEnd%); 'get the mean and SD of
baseline calcium activity
325     BaseSD[i%] := ChanMeasure(i%+1, 12, BaseStart, BaseEnd%);
326     Thresh[i%] := BaseMean[i%] + BaseSD[i%] * Mult; 'calculate threshold for calcium
event detection
327 next
328
329
330 end
331
332
333 ''''''''''''''Create event chan, mark suprathreshold events detected in each source
wavechan
334 proc ThreshMark(mt%, Thresh[], NumWav%, MinDur%, &EventChan%[], BaseEnd%);
335
336 var mc%, cc%; 'memory channel index, counter of suprathreshold events per channel
337 var i%, ii%; 'loop indices
338
339 for i%:=0 to NumWav%-1 do 'for all cells
340

```

```

341     cc%:=0; 'reset suprathreshold counter
342     mc%:=MemChan(4, 0); 'create temporary chan for storing of suprathreshold event
periods (level)
343
344     'Scan all data points after baseline, see where thershold is crossed, then see if the
minimal duration rule is fulfilled;
345     'Mark in level type channel with high level indicating event duration:
346
347     for ii%:=BaseEnd% to mt%-1 do 'go over file, from end of baseline until the end, look
for continuous suprathreshold events
348         if ChanMeasure(i%+1, 2, ii%-0.25, ii%+0.25) > Thresh[i%] then 'suprathreshold
point?
349             cc% += 1; 'add to counter
350
351             if cc% = MinDur% then 'duration rule fulfilled?
352                 MemSetItem(mc%, 0, ii%-MinDur%+1); 'mark level transition from low to high
353             endif
354
355             else 'value at this timepoint is below threshold
356                 if cc% >= MinDur% then 'were previous timepoints within a detected event?
357                     MemSetItem(mc%, 0, ii%); 'mark level transition from high to low
358                 endif
359
360             cc%:=0; 'reset suprathreshold counter
361         endif
362     next
363
364     EventChan%[i%] := ChanSave(mc%, i%+NumWav%+1); 'save to permanent channel and store
identifier
365     ChanShow(EventChan%[i%]); ChanTitle$(EventChan%[i%], "Events " + Str$(i%+1)); 'name
event chan to identify the waveform chan it refers to
366     ChanDelete(mc%); 'delete memory chan
367
368 next
369
370 end
371
372
373
374 '::::::::::::::::::::::::::'For each calcium event, find peak amplitude and area over threshold
375 proc EventParams(&RiseTime[[]], &FallTime[[]], &EventDur[[]], &Peak[[]], &AOT[[]],
NumWav%, EventChan%, mt%, MinDur%, Thresh[], BaseEnd%);
376
377 var et[NumWav%][mt% - BaseEnd%]; 'temporary array for eventtimes
378
379 var i%,ii%; 'loop indices
380 var c%,j%,k%; 'counters
381
382 for i%:=0 to NumWav%-1 do 'for all channels
383     c% := Count(EventChan%[i%], BaseEnd%, MaxTime()); 'count rising/falling level events
per chan
384
385     if c% mod 2 <> 0 then 'missing falling event (over thershold until the end of recording)
386         c% += 1;
387     endif
388
389     j% := 0; 'reset rise/fall counters
390     k% := 0;
391
392     if c% > 0 then 'are there any suprathreshold events in this channel?
393         ChanData(EventChan%[i%], et[i%][], BaseEnd%, MaxTime()); 'find rising/falling times
394
395         for ii%:=0 to c%-1 do 'go over all rising/falling, group to rising and falling
396
397             if (ii% mod 2) = 0 then 'rising events
398
399                 RiseTime[i%][j%] := et[i%][ii%];
400                 j%+=1; 'add to rising events counter
401
402             else 'falling events
403
404                 FallTime[i%][k%] := NextTime(EventChan%[i%], et[i%][ii%-1]);
405
406                 if FallTime[i%][k%] < 0 then 'if event doesn't go down until end of
recording (returns negative value)
407                     FallTime[i%][k%] := MaxTime(); 'consider end of recording as end of

```

```

event
    endif
408
409
410         k%+=1; 'add to falling events counter
411     endif
412
413     next
414
415     for ii%:=0 to c%/2-1 do 'measure parameters between each rise and fall
416         EventDur[i%][ii%] := FallTime[i%][ii%] - RiseTime[i%][ii%]; 'event duration
417         Peak[i%][ii%] := ChanMeasure(i%+1, 8, RiseTime[i%][ii%], FallTime[i%][ii%]);
'peak amplitude of event
418         AOT[i%][ii%] := ChanMeasure(i%+1, 4, RiseTime[i%][ii%], FallTime[i%][ii%]);
'event area (over threshold)
419         AOT[i%][ii%] := AOT[i%][ii%] - (Thresh[i%] * EventDur[i%][ii%]);
420     next
421
422     endif
423
424 next
425
426
427 end
428
429
430 '::::::::::::::::::::::::::'print txt file with the drift-corrected values per chan per s, same
format as original txt files
431 proc PrintResults_dFF(mt%, vh%, NumWav%, FileEnd1$, &t1%, CloseWins%, OGFN$);
432
433 FrontView(vh%); 'go to data file
434
435 var ValperSec[NumWav%][mt%]; 'matrix for storing [filtered] dF/F values per cell per sec
436
437 var i%,ii%; 'loop indices
438
439 for i%:=0 to NumWav%-1 do 'for each cell
440     for ii%:=0 to mt%-1 do 'for each time point
441         ValperSec[i%][ii%]:=ChanMeasure(i%+1, 2, ii%, ii%+0.99); 'corrected value
442     next
443 next
444
445 var vps$; 'for filling rows w text
446 for ii%:=0 to mt%-1 do 'fill rows one by one
447     for i%:=0 to NumWav%-1 do 'go column by column
448         if i% = NumWav%-1 then 'if end of row, go to next row; otherwise move to next
column
449             vps$ += (Str$(ValperSec[i%][ii%]) + "\n");
450         else
451             vps$ += (Str$(ValperSec[i%][ii%]) + "\t");
452         endif
453     next
454 next
455
456 FileNew(1, 1); 'new txt file and view handle
457 t1%:=FrontView();
458
459 Print(vps$); 'print all columns and rows into txt file
460
461 var NFN$; 'name and save new file
462 NFN$:= OGFN$ + "_" + FileEnd1$;
463 FileSaveAs (NFN$, 1);
464
465 if CloseWins% = 1 then
466     FileClose(0, 0);
467 endif
468
469
470 end
471
472
473
474 '::::::::::::::::::::::::::'create csv file with peak, area, threshold, and summary for all cells in
this recording
475 proc PrintResults_Summary(vh%, HPFreq, Mult, NumWav%, MaxBins%, EventDur[[]], Peak[[]],
AOT[[]], Thresh[], MinDur%, FileEnd2$, DCTC, CloseWins%,
476 OGFN$, FT$);

```



```

477
478
479 var ef%; 'indicator for if there are no suprathreshold events ("empty file")
480 var RealMaxBins%; 'number of most detected events in single channel - for determining
    number of rows
481
482 FrontView(vh%); 'go to data file
483
484 var C$, E$, P$, A$, T$, I$; 'for storing strings related to each parameter
485
486 C$:="Channel (cell) No.:\n";
487 E$:="Event durations (s):\n\n";
488 P$:="Peaks (absolute values):\n\n";
489 A$:="Areas (over threshold):\n\n";
490 T$:="Threshold per channel:\n";
491
492 if DCTC > 0 then 'if DC remove process was used, indicate this in summary
493
494     I$:="Analysis info:\nIIR drift correction frequency = "+Str$(HPFreq)+" ("&FT$&"),
    followed by DC remove ("&Str$(DCTC)+
495     " s)\nThreshold = mean + " &Str$(Mult)+" x SD of baseline\nMinimal event duration =
    "&Str$(MinDur%)+&" s";
496
497 else
498     I$:="Analysis info:\nIIR drift correction frequency = "+Str$(HPFreq)+"
    ("&FT$&")\nThreshold = mean + " &Str$(Mult)+
499     " x SD of baseline\nMinimal event duration = "&Str$(MinDur%)+&" s";
500
501 endif
502
503
504 var i%,ii%,c%; 'loop indices, counter
505
506 for i%:=0 to NumWav%-1 do 'go over all channels
507     C$ += Str$(i%+1) + ", ";
508     c% := 0; 'initialize event counter
509     for ii%:=0 to MaxBins%-1 do 'find maximum number of events across all channels
510         if Peak[i%][ii%] > 0 then
511             c% += 1; 'add event to count
512         endif
513     next
514     if c% > RealMaxBins% then 'new highest event count
515         RealMaxBins% := c%;
516     endif
517
518
519 next
520
521 docase
522 case RealMaxBins% = 0 then ef% := 1 'no calucium events detected
523 case RealMaxBins% < 3 then RealMaxBins% := 3; 'minimum of 3 lines space
524 endcase
525
526
527 for ii%:=0 to RealMaxBins%-1 do 'go over all peak values per row
528     for i%:=0 to NumWav%-1 do 'go over all channels
529         if i%=NumWav%-1 then 'if last column, move to next row
530             if EventDur[i%][ii%]>0 then 'event detected? if not, skip to next
531                 E$+=Str$(EventDur[i%][ii%])+"\n";
532                 P$+=Str$(Peak[i%][ii%])+"\n";
533                 A$+=Str$(AOT[i%][ii%])+"\n";
534             else
535                 E$+="\n";
536                 P$+="\n";
537                 A$+="\n";
538             endif
539
540
541         else
542             if EventDur[i%][ii%]>0 then 'event detected? if not, skip to next column
543                 E$+=Str$(EventDur[i%][ii%])+", ";
544                 P$+=Str$(Peak[i%][ii%])+", ";
545                 A$+=Str$(AOT[i%][ii%])+", ";
546             else
547                 E$+=", ";
548                 P$+=", ";

```



```

549         A$+=",";
550     endif
551 endif
552
553 next
554 next
555
556 for i%:=0 to NumWav%-1 do 'threshold value per channel
557     T$+=Str$(Thresh[i%])+",";
558 next
559
560 FileNew(1, 1); 'new text file to print into (to be saved as CSV)
561
562 if ef% = 1 then
563     Print("*****No suprathreshold events detected in file!*****\n\n\n", T$, "\n"*3,
564     I$);
565 else
566     Print(C$, "\n"*3, E$, "\n"*3, P$, "\n"*3, A$, "\n"*3, T$, "\n"*3, I$);
567 endif
568
569 var NFN$;
570 NFN$ := OGFN$ + "_" + FileEnd2$; 'save and close csv file
571 FileSaveAs(NFN$,1); FileClose(0, 0);
572
573
574 end

```